smartLAB[®]

fit W

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smartLAB[®]

fit W

Body Analyser Scale with Bluetooth Smart and ANT+

User Manual

Please read the user manual carefully before the first use.

Current information about this device, visit smartLAB on

www.support.hmm.info

HMM Diagnostics GmbH
Friedrichstraße 89
D-69221 Dossenheim, Germany
www.hmm.info

Content

I. Introduction	25
Principle of measurement	25
II. Your smartLAB[®]fit W	26
Display & functions	26
Specifications	27
Contents	27
General Information	27
Security information	29
III. Setup & Operating functions	30
Inserting and changing the batteries	30
Beginning of operation / „quick weighing “	30
Setting up of the user data	31
Function setup	33
Measuring	34
Measuring only weight	35
Wireless data transfer	35
Data transfer via user profile	36
Data transfer to other wireless data transmitting devices	37
Optional: Weighing on a carpet	37
Interpreting the results	38
IV. Miscellaneous	41
Maintenance	41
Regularity and Safety Notice	42
Warranty	43
Affichage & fonction	47

I. Introduction

Thank you for having decided in favor of our **smartLAB[®]fit W** Body Analyser Scale. This product enables you to measure your body weight as well as content of adipose, content of water, content of muscle and bone mass. Please read the instructions carefully before you first use the scale and keep the instructions for a future use.

Principle of measurement

Your **smartLAB[®]fit W** Body Analyser Scale is based on the principle of the bioelectric impedance analysis (BIA). A weak and non-perceptible and, therefore, not dangerous electric signal will be conducted through the whole body. The scale identifies the electric resistance (impedance). Individual data such as age, body height, gender or degree of activity allow to identify the adipose content and other parameters.


Muscle tissue and water have a good electric conductivity and, therefore, a lower resistance than bone- and adipose tissue because adipose cells and bones hardly conduct current.

Please note that the values of our **smartLAB[®]fit W** Body Analyser Scale are only an approach to the real medical values of the body. Only a specialist who takes medical methods into account (e. g. computer tomography) can identify the exact parameters like content of adipose, content of water, content of muscles and bone mass.


The **smartLAB[®]fit W** Body Analyser Scale is an ANT+ and *Bluetooth[®]* Smart device. It supports the ANT+ device profile and is compatible with receiving devices (such as sport watches, PC and Smartphones) that support the ANT+ weight scale device profile, e.g. the Samsung S Health App supports the connection of ANT+ Scales. With the integrated *Bluetooth[®]* Smart the **smartLAB[®]fit W** Body Analyser Scale can also sent data via Bluetooth Smart to a corresponding receiver. You can then transmitt your values wireless to these receiving devices. The **smartLAB[®]fit W** Body Analyser Scale will choose the correct data transfer technology for you.

II. Your smartLAB[®]fit W

Display & functions

The settings can be adjusted by pressing on the side of the scale.

Specifications

1. *Device: smartLAB[®]fit W*
2. *Max. capacity: 150 kg / 330 lb*
3. *Units: kg / lb (in 100 g / 0,2 lb steps)*
4. *Buttons: Set, Go up, Go down*
5. *Activity indices: 4*
6. *Adipose: in % for 0,1% steps*
7. *Content of water of the body: in % for 0,1 % steps*
8. *Muscle mass: in % for 0,1% steps*
9. *Basal metabolic syndrome: in kcal with 1 kcal per step*
10. *Active metabolic syndrome: in kcal with 1 kcal per step*
11. *User profiles: 8*
12. *Battery: 4 x AAA batteries (1,5 V)*
13. *Size (L x W x H): 325x325x19*
14. *Weight: 1950 g*
15. *Operating temperature: 10°C-40°C*
16. *Storage temperature: -20°C - 60°C*
17. *Symbol when overburdened: „----“*
18. *Symbol for low batteries: „Lo“*
19. *Communication/Data transfer: ANT+ and BLE wireless*

Contents

1 smartLAB[®]fit W Body Analyser Scale

1 User manual

4 1,5V „AAA“ batteries

4 Carpet extension feet

General Information

- Please place the scale on an even, hard floor. The floor is the basis for a correct measurement.
- If possible please use the scale at the same time every day (we recommend in the morning), after having used the toilet, with empty stomach and without clothes in order to be able to compare the results.
- Please wait about 15 minutes after getting up before you measure your body weight so that the water of your body can disperse.
- The adipose can only be detected bare-footed. You can moisten the soles of your feet. Completely dry feet can result an incorrect measurement because they do not have enough conductivity.
- Please do not move during the measurement.
- Please wait several hours before measuring after unaccustomed physical exercise.

- The following persons or groups of persons will not get differing or not plausible results concerning their adipose and the other values:
 - children under 10 years
 - competitive athletes & bodybuilders
 - pregnant woman
 - people having fever, dialysis patients, people with edema symptoms or osteoporosis
 - persons taking cardiovascular medicine
 - persons taking vasodilatory or vasoconstrictive medicine
 - persons with considerable anatomic differences of their legs that can influence the height of the body (different lengths of the legs).
- The device is for your own use only and not for medical or commercial use.
- Please note that technical related measuring tolerances are possible since the scale is not especially made for professional medical use.
- The **smartLAB[®] fit W** Body Analyser Scale corresponds to the EC directive 2004/108.

Please feel free to contact the manufacturer or the customer service for more information.

Security information

The **smartLAB[®]fit W** Body Analyser Scale must not be used by persons with medical implants (e. g. cardiac pacemaker). Their function could be affected.

Do not use during pregnancy. The amniotic liquor inside the body can influence the measurements.

Please do not use the **smartLAB[®]fit W** Body Analyser Scale with wet feet – slip hazard!

Do not swallow the batteries – danger of life. Please store batteries and Body Analyser Scale at a place children cannot reach. In case a battery has been swallowed, please consult a physician immediately.

Keep children away from the packaging materials - danger of asphyxiation!

Do not recharge the batteries. Do not take them to pieces. Do not throw them into a fire. Do not short-circuit them.

III. Setup & Operating functions

Inserting and changing the batteries

- Please open the cover of the battery case at the back of the device.
- Insert 4 AAA batteries (1,5 V). Please pay attention to the polarity.
- Close the cover of the battery case.

Important: Please make sure that your scale is close to a gateway/App when inserting and changing the batteries, so the gateway/App can set the scale and send important information like date, time and profile information. The „W“-Symbol on the scale indicates when the scale is trying to communicate with the gateway.

Note: In case the scale does not operate please remove the batteries and insert them again.

Your **smartLAB[®]fit W** Body Analyser Scale has a symbol indicating that the batteries need to be changed. In case the batteries are almost empty, the symbol “Lo” can be seen on the display. You should change the batteries then.

Please remove the batteries in case you are not going to use the scale for one month or more in order to avoid an automatic discharge of the batteries.

The Body Analyser Scale as well as the batteries must be disposed of according to the local regulations.

In order to guarantee a correct measurement you should calibrate the scale always after any location change. Just tip on the scale and wait until it shows “0.0” on the scale.

Beginning of operation / „quick weighing “

- Use the **smartLAB[®]fit W** Body Analyser Scale bare-footed. Do not move so you can stand on both legs. The scale starts measuring immediately

and finally shows your body weight on the display. Now you can leave the tread of the scale. The scale will deactivate automatically after several seconds.

Setting up of the user data

In order to measure the adipose content and other parameter you first have to insert your individual user data. You can store up to 8 user profiles and recall the measured data.

The **smartLAB[®]fit W** Body Analyser Scale has a special function for the automatic transfer of user profile information and function setup information. It is required that the receiving devices supports these functions. This automatic setup function is supported by Apps and Gateways by HMM.

With the automatic setting of the scale there is no need to setup the scale manually. When you use the scale for the first time or you change the batteries, the scale will check for an adequate receiving device. Profile and setup information is also transferred after each data transmission to the receiving device.

If a receiving device is in reach of the scale, it will send all the setup data to the scale. This function has the benefit of also sending the batterie status to the receiving device.

The automatic setting can be disabled with the function setup (Heart Beat) (please refer to the section function setup).


If you want to set up the scale manually please press "SET" (on the side of your scale) to create a new user profile.

Note: The weight scale can be set by the touch function, which can be found on the right top side of the scale. The touch function represents the buttons, which are abreast of the Set writing and the arrows at the side of the scale.

User profile:

The display now shows several illuminated preadjustments. The symbol "P1" appears on the display (bottom, left) for user profile 1 (see pict. 1).

You can confirm this choice by pressing “SET” or you can create a new user profile by increasing the number (arrow keys on the side of the scale).

Pic. 1

Gender:

Press „SET“ to confirm your choice and then chose your gender. The display shows the symbol for a man. If you want to change the gender please use one of the arrow keys until you can see the symbol for a woman. Press “SET” to confirm your choice.

Body height:


Please insert your body height by using the arrow keys. You can insert body heights from 100 cm up to 250 cm. When you continue to press the buttons the indicated values go faster. Please confirm your value by pressing “SET”.

Age:

The display shows the value for your age. Please use again the arrow keys to insert your age and confirm by pressing “SET”.

Activity index:

Then the symbol for the choice of the activity index appears in the left display segment (see picture 2).

Pic. 2

The value defines a certain degree of activity that is important for a medium-term or a long-term consideration. You can choose between 4 activity degrees for each profile by using the arrow keys. Please confirm your choice and press “SET”.

Defining the different degrees of activity:

- Degree 1: No physical exercise
- Degree 2: Only little exercise such as promenades, working in the garden or gymnastics.
- Degree 3: medium physical exercise (30 minutes physical exercise from 2 up to 4 times per week minimum)
- Degree 4: high physical exercise (30 minutes physical exercise from 4 up to 6 times per week minimum)


Please note that the fat calculation is depending on the activity index. The higher the index the lower the fat calculation will be. This has physical reasons.

Function setup


After the input of the activity index you will enter the function setup. The function setup can be set to ON or OFF.

Heart Beat : This is the automatic setting function. Switch through the arrows to ON or OFF and then confirm with „SET“.


Pic. 3

Beep 1 : A long beep will be heard if the measurement is finished. Switch through the arrows to ON or OFF and then confirm with „SET“.


Pic. 4

Beep 2 : 3 short beeps for transferring only weight or 5 short beeps for transfer of all data. This indicates that the transmission was successful. Switch through the arrows to ON or OFF and then confirm with „SET“.

Pic. 5

Beep 3 : 2 long beeps will be heard if the transmission was NOT successful. Switch through the arrows to ON or OFF and then confirm with „SET“.

Pic. 6

After having finished the setting up of your data the symbol “0,0” appears on the display. The **smartLAB[®]fit W** Body Analyser Scale is ready to start the measurement.

Note: Profile and function setup will only be saved on the scale (when setting manually) when going through the settings until the display shows „0,0“. If not going through the setup until the end, the entered data information will not be saved on the scale.

Note: The beep tones are turned off by factory settings. These can either be turned on manually or automatically by a HMM App or Gateway.

Measuring

After having determined all parameters and user profiles body weight, adipose and the other values can be determined.

- Please press „SET“ and choose your user profile by using the arrow keys. Do NOT confirm by pressing “SET” since you would have to insert and confirm again all parameters such as gender, body weight etc.
- The measurement starts automatically after a few seconds. First of all your **smartLAB[®]fit W** Body Analyser Scale detects your body weight (left display segment). After that the measurement of the other parameters will start. You will see three circles on the display segments concerned (please refer to chapter “display & functions” in this user manual). The scale is detecting the other parameters – content of water, content of muscles and bone mass (see picture 7).

Pic. 7

- Finally all values are displayed (see picture 8). After a few seconds the indication of the body weight changes into the indication of the calorie consumption (see picture 9), then into the indication of the activity (in kcal). In addition the activity index is indicated by the symbol "activity index" (see picture 10). At the end your **smartLAB[®]fit W** Body Analyser Scale displays again the parameters for your user profile which are the basic data for your measured values.

The **smartLAB[®]fit W** Body Analyser Scale recognize automatically through your weight, which person is currently standing on the scale and uses his user profile. This automatic function requires that you once select the desired user and weight yourself barefoot after that. Future measurements will be recognized automatically.

Measuring only weight

When you step on the scale not bare-footed the scale will measure the weight only.

Wireless data transfer

Your **smartLAB[®]fit W** Body Analyser Scale has an integrated ANT+ and Bluetooth Smart module which enables you to transfer measured data to your computer, laptop, PC, Smartphone or sport watch. There you can save, analyze and process the values

Please check our website www.hmm.info for a list of Apps that work with the scale under „Apps for smartLAB product“.

When the scale has been successfully paired with a receiving device, the measured data will automatically be sent to the corresponding device such as Smartphone or PC via ANT+ or Bluetooth Smart. The **smartLAB[®]fit W** Body Analyser Scale will automatically select the correct data transfer technology.

When using the HMM Apps (such as fitmeFit move App or hLine App) or Gateways, the profile and function settings can be automatically transmitted to the scale. This happens during a data transfer or after a battery change/ inserting the batteries. Please find more information in the chapter „Setting up the user data“.

Please find more information on the connection possibilities and other features of the scale under www.support.hmm.info.

Note: HMM Apps are the fitmeFit move App and hLine App, both available for Android and iOS.

What is ANT?

ANT is a wireless data transfer standard. In order to use it, both transmission- and receiving device must be equipped with an ANT module.

What is BLE?

BLE is a wireless data transfer standard and is part of *Bluetooth[®] Smart*. In order to use it, both transmission- and receiving device must be equipped with an BLE module.

Data transfer via user profile

Set up your user profile according to the instructions in the “Measuring” chapter of this manual. Your **smartLAB[®]fit W** Body Analyser Scale is preset to transfer the data automatically. Please choose your user profile as it is described in the chapter “Measuring” and start the measurement. Your **smartLAB[®]fit W** Body Analyser Scale will transfer the measured data to the receiving device automatically after the measurement.

It is important that the correct user is used when transferring data. If the scale has not been used before, select your user

profile by pressing the “Set” button and the desired user through the up and down arrows. After weighing the scale saves this weight to that user.

The **smartLAB[®]fit W** Body Analyser Scale recognizes automatically through your weight, which person is currently standing on the scale and uses his user profile. This automatic function requires that you once select the desired user and weigh yourself after that. Future measurements will be recognized automatically.

Data transfer to other wireless data transmitting devices

You can connect your **smartLAB[®]fit W** with a partner device which fulfills the standard protocol of ANT or Bluetooth Smart wireless data transfer. You can transfer measured weight values automatically and wireless to the device. Partner devices could be laptops, smartphones, tablets or any other ANT or Bluetooth Smart-capable device. The communication will be enabled by a gateway* which can be installed on PC, Smartphone or Tablet. Make sure that the correct profile is used when weighing.

Optional: Weighing on a carpet

You can use the **smartLAB[®]fit W** Body Analyser Scale with the extension feet on a carpet. For this please follow these steps:

1. Remove the non-slip pads from the scales feet

2. Fit the extension feet firmly onto the scales feet. Now you can use the scale on a carpet.

Interpreting the results

The following guiding values help you to interpret the measured values of your **smartLAB®fit W** Body Analyser Scale. These values are only an approach you should eventually discuss and interpret with your physician.

Adipose content in %:

The following values are guiding values concerning your adipose content. Please contact your physician for further information.

Age	Woman				Man			
	undernourished	healthy / normal	overweight	too much adipose	undernourished	healthy / normal	overweight	too much adipose
10 - 12	< 12	12 - 22	22 - 30	> 30	< 8	8 - 18	18 - 24	> 24
13 - 18	< 15	15 - 25	25 - 33	> 33	< 8	8 - 18	18 - 24	> 24
19 - 30	< 20	20 - 29	29 - 36	> 36	< 8	8 - 18	18 - 24	> 24
31 - 40	< 22	22 - 31	31 - 38	> 38	< 11	11 - 20	20 - 26	> 26
41 - 50	< 24	24 - 33	33 - 40	> 40	< 13	13 - 22	22 - 28	> 28
51 - 60	< 26	26 - 35	35 - 42	> 42	< 15	15 - 24	24 - 30	> 30
60+	< 28	28 - 37	37 - 47	> 47	< 17	17 - 25	26 - 34	> 34

Table 1

Note: Athletes often have lower values than those of the table above. The values can be underneath the guiding values depending on the sport, the intensity of exercise and the physical condition. But please note that extremely low values can be an indication for health hazard.

Content of water in your body in %:

The content of water inside your body is usually within the following ranges (only the normal range is indicated):

Age	Woman	Man
10 - 12	> 60	> 64
13 - 18	> 58.5	> 63.5
19 - 30	> 56	> 62.5

Age	Woman	Man
31 - 40	> 53	> 61
41 - 50	> 52	> 60
51 - 60	> 51	> 59
60+	> 50	> 58

Table 2

Note: Adipose has relatively few water. For this reason the content of water in case of a high content of adipose can be below the guiding values above. Regarding endurance the values of others can be higher than the guiding values because of the low adipose and the high content of muscles.

Content of muscles in %:

The content of the muscles usually stays within the following ranges (only the normal range is indicated):

Age	Woman	Man
10 - 99	> 34	> 40

Table 3

Bone mass in kg:

Our bones as well as our whole body are subject to the ageing process and all its consequences. During the childhood the bone mass is increasing fast and reaches its maximum at the age of 30-40. Then the bone mass is decreasing. You can influence the ageing process by a healthy nutrition (esp. calcium and vitamin D) and by a continuous physical exercise. In addition to that you can sustain your bones by a specific muscle formation.

Bone mass and bone density are different factors. Please do not confuse them. The latter can only be detected by a medical examination (e. g. computer tomography, ultra sound). For this reason the scale cannot be used to identify changes or the bones or their durability.

The bone mass can hardly be influenced but is changing according to factors like weight, body height, age or gender.

	Woman			Man		
Body weight	< 45 kg	45 - 60 kg	> 60 kg	< 60 kg	60 - 75 kg	> 75 kg
Bone mass	1,8 kg	2,2 kg	2,5 kg	2,5 kg	2,9 kg	3,2 kg

Table 4

Basal Metabolic Rate (BMR):

BMR means the amount of energy the body needs for its basal vital functions when the body is in the state of absolute calmness (e. g. when you stay in bed for 24 hours). This value mainly depends on factors like weight, body size and age. This value is displayed on your **smartLAB®fit W** Body Analyser Scale. The unit is kcal/Tag (kcal/day). The basis for the calculation is the science-based Harris-Benedict-Formula.

This is the amount of energy your body needs at least in the form of nutrition. If your body gets less energy for a longer time this can affect your health negatively.

Active Metabolic Rate (AMR):

The active metabolic rate means the amount of energy the body needs daily in the state of activity. The energy consumption increases with an increasing physical exercise. Your **smartLAB®fit W** Body Analyser Scale has 4 different activity indices (refer to chapter "Setting up of the user data") to determine the appropriate degree of activity.

In order to maintain the current weight the consumed energy must be replaced in form of food and drink. In case this energy will not be replaced the body uses the stored fat – the body loses weight. In case the body gets more energy than the required AMR the body is not able to use this energy – it will be stored inside the body and you gain weight.

Time as a factor:

Please note that only a long-term result counts. Short-term changes of your weight within a few days are absolutely normal and are often due to changes of the content of water. Only middle-term or long-term changes are related to the content of your adipose or your muscles. In case your weight decreases for a short time and the content of the adipose of your body increases or stays the same you only lost water (e.g. after a physical exercise or after having visited a sauna). If your weight increases in the medium term and your content of adipose decreases or stays the same you might have build up more muscle mass.

Note: Adipose, content of water or muscle mass may not be added because the muscle tissue also contains components of water.

IV. Miscellaneous

Maintenance

- Please clean your device regularly by a damp cloth. Do not use caustic cleaning solutions. Do not put the device under fluent water.
- Protect the device from impacts, humidity, dust, chemicals or fluctuations of temperature. Do not store near heat sources such as ovens or heating installations.
- NEVER take the scale to pieces by yourself. In this case the warranty will not be valid anymore. Only the manufacturer or the authorized trader is allowed to repair the scale. Please check the batteries before any complaint and change them if necessary.
- The **smartLAB[®]fit W** Body Analyser Scale is in accordance with the EC directive 2004/108. Please feel free to contact the manufacturer or the customer service for further questions regarding the applications of the scale.

Regularity and Safety Notice

Declaration of Conformity

This product is compliant with the EC Directives: -1999/5/EC

Safety Statement

This product has been verified to conform to the safety certification in accordance with the specification of EN Standards: EN 60950-1:2006 + A11:2009 + A1:2010 + A12:2011

WEEE note

The WEEE (Waste Electrical and Electronic Equipment) Directive, that came into effect as European law on 13th February 2003, led to a major change in the disposal of electrical equipment.

The primarily purpose of this directive is to prevent electronic waste (WEEE), while encouraging for reuse, recycling and other forms of reprocessing to reduce waste.

The WEEE logo on the product and the packaging indicates that it is not allowed to dispose of the product in the ordinary household waste. It is your responsibility to deliver all the disused electrical and electronic equipment to the respective collection points. A separate collection and proper recycling of electronic waste helps dealing economically with natural resources. Furthermore, the recycling of electrical waste is a contribution to the conservation of our environment and thus the health of all people. More information concerning the disposal of electrical and electronic devices, reprocessing and the collection points you can get from the local authorities, waste disposal companies, from retailers and manufacturers of the device.

RoHS compliance

This product complies with Directive 2002/95/EC of the European Parliament and the Council of 27th January 2003 regarding the limited use of dangerous substances in electrical and electronic equipment (RoHS) and its variations.

Warranty

HMM Diagnostics GmbH products need to fulfill high quality requirements. Because of this reason, HMM Diagnostics GmbH gives a 2-year warranty by purchasing this **smartLAB**[®] product. You can even extend the warranty from 3 to 5 years without extra pay when you register your product. Please register online under following website: www.hmm.info/en/registration.

Wear parts, batteries etc. are excluded from warranty

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*The Bluetooth[®] word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by HMM Diagnostics is under license. Other trademarks and trade names are those of their respective owners.