smartLAB[®] pressure

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Upper Arm Blood Pressure Monitor

User Manual

Please read the user manual carefully before the first use.

Current information about this device, visit smartLAB on www.support.hmm.info

HMM Diagnostics GmbH
D-69221 Dossenheim, Germany
www.hmm.info

Content

I. Introduction

Product description	36
Normal variations of your blood pressure	36
Important Information	37

II. Your smartLAB[®] *pressure*

Device	40
Specifications	42
Contents	42
General Information:	43

III. Setup & operating functions

Making the device ready for operation	44
Setting time, date and average	44
Select the User	46
Connecting the cuff to the monitor	47
Applying the cuff	47
Body posture during measurement	48
Measurement	49
Displaying stored values	50
Deleting measured values	52
Lock the buttons	53

IV. Miscellaneous

Maintenance	54
Explanation of symbols	54
Troubleshooting	55
Assessing high blood pressure for adults	56
DPA technology for a more accurate measurement	60
Measure-during-inflation	60
Regularity and Safety Notice	61
Warranty	63

36

36

37

40

40

42

42

43

44

44

44

46

47

47

48

49

50

52

53

54

54

54

55

56

60

60

61

63

I. Introduction

Thank you for having decided for our **smartLAB® pressure** blood pressure monitor for the upper arm. Our device enables you to control your blood pressure and your pulse rate very easily. Furthermore, you can store your values by an internal measured data storage. The classification by colors of the measured values on the side of your device enables you to interpret your values immediately.

Please read carefully the instructions before you use the **smartLAB® pressure** blood pressure monitor for the first time.

Product description

Blood pressure and heart rate can be identified automatically and without any invasive method by an integrated pressure sensor technology and by an oscillometric measurement method.

Before every measurement, the unit establishes a “zero pressure” equivalent to the air pressure. Then it starts inflating the arm cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate.

The device also compares the longest and the shortest time intervals of detected pulse waves to mean time interval then calculates standard deviation. The device will display a warning signal with the reading to indicate the detection of irregular heartbeat when the difference of the time intervals is over 25%.

Blood pressure, pulse rate and (possible) irregular heartbeats can be seen on the LCD display. You can store the last 60 results (2 user profiles) with date and time. In addition to that you can see the average of the last 3 measurements.

Normal variations of your blood pressure

Many factors influence your blood pressure such as physical exercises, excitement, stress, nutrition, drinking, smoking and other activities (including measuring your blood pressure). For this reason you almost never get the same results.

Your blood pressure is subject to permanent variations – day and night. The highest values are usually in the day, the lowest usually at midnight. Usually your blood pressure starts

increasing after 3 a.m. and reaches its highest value in the daytime when you are awake and active.

For this reason we recommend to measure your blood pressure every day at the same time under the same circumstances. Please relax 3 to 5 minutes between the measurements so that your blood circulation can start again in your arm. It is very rare that you get identical measured values when you measure your blood pressure several times subsequently.

Important Information

CAUTION


- The device is not suitable for measuring the blood pressure of children. Ask your doctor before using it on older children.
- The device is not suitable for use on pregnant women, patients with implanted, electrical devices, patients with pre-eclampsia, premature ventricular beats, atrial fibrillation, peripheral, arterial disease and patients undergoing intravascular therapy or arterio-venous shunt or people who received a mastectomy. Please consult your doctor prior to using the unit if you suffer from illnesses.
- Do not take any therapeutic measures on the basis of a self measurement. Never alter the dose of a medicine prescribed by a doctor. Consult your doctor if you have any question about your blood pressure.
- Please keep the unit out of reach of infants, children or pets, since inhalation or swallowing of small parts is dangerous or even fatal. This device is intended only for adult use at home.
- This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than arm or for functions other than obtaining a blood pressure measurement.
- If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the START/STOP button to release the air immediately from the cuff. Loosen the cuff and remove it from your arm.
- On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, open the cuff immediately. Prolonged high pressure (cuff pressure >300mmHg or constant pressure >15mmHg for more than 3 minutes) applied to the arm may lead to an ecchymosis.

- Too frequent and consecutive measurements could cause disturbances in blood circulation and injuries.
- Do not wrap the cuff on the same arm which other monitoring ME EQUIPMENT is applied simultaneously, because this could cause temporary loss of function of those simultaneously-used monitoring ME EQUIPMENT.
- Don't kink the connection tube during use, otherwise, the cuff pressure may continuously increase which can prevent blood flow and result in harmful injury to the PATIENT.
- The equipment is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- This unit is not suitable for continuous monitoring during medical emergencies or operations.
- This device cannot be used with HF surgical equipment at the same time.
- This device is not intended for patient transports outside a healthcare facility.
- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- The operator shall not touch output of batteries/adaptor and the patient simultaneously.
- Manufacturer will make available on request circuit diagrams, component parts list etc.
- The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010.
- It will not cause any potential sensitization or irritation reaction. Never apply the cuff over hurt skin.
- Do not use the unit in case of existing polyester resp. synthetic allergies.
- Be careful to strangulation due to cables and hoses, particularly due to excessive length.
- Do not connect the air hose to other medical equipment, as this could cause air to be pumped into intravascular systems or high pressure, what could lead to dangerous injuries.


- Before use, make sure the device functions safely and is in proper working condition.
- Please use the device under the environment which was provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.
- Please use ACCESSORIES and detachable partes specified/ authorised by MANUFACTURER. Otherwise, it may cause damage to the unit or danger to the user/patients.
- Please dispose of ACCESSORIES, detachable parts, and the ME EQUIPMENT according to the local guidelines.
- Please do not attempt to repair the unit yourself in the event of malfunctions. Only have repairs carried out by authorized service centers.
- Please report to Manufacturer if any unexpected operation or events occur.
- The device doesn't need to be calibrated in two years of reliable service.
- Please use the soft cloth to clean the whole unit. Don't use any abrasive or volatile cleaners.

II. Your smartLAB® pressure


Device

Back side

LCD display

SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic blood pressure	High blood pressure
DIA	Diastolic blood pressure	Low blood pressure
PUL /min	Pulse display	Pulse in beats per minute
AVG	Average value	The average value of the latest three records
Q	Memory	Indicate it is in the memory mode
kPa	kPa	Measurement Unit of the blood pressure (1kPa=7.5mmHg)
mmHg	mmHg	Measurement Unit of the blood pressure (1mmHg=0.133kPa)
Lo + 	Low battery	Batteries are low and need to be replaced
	Irregular heartbeat	Blood pressure monitor is detecting an irregular heartbeat during measurement.
 	Blood pressure level indicator	Indicate the blood pressure level
88/88	Time	Year/Month/Day, Hour : Minute

SYMBOL	DESCRIPTION	EXPLANATION
	Motion indicator	Motion may result in an inaccurate measurement.
	Heartbeat	Blood pressure monitor is detecting a heartbeat during measurement.
	User 1	Start measurement and save the measuring results for User 1
	User 2	Start measurement and save the measuring results for User 2

Specifications

1. Device type: **smartLAB[®]pressure**
2. Classification: class II, type B
3. Size of device: 107 mm x 103 mm x 118 mm
4. Size of cuff: 420 mm x 200 mm
5. Weight: 250 g (without batteries and Cuff)
6. Measurement method: oscillometric, automatic inflation and measurement
7. Storage capacity: 60 values with date and time (2 user profiles)
8. Batteries: DC 6V or 4 „AAA“ batteries or alternative AC adapter
9. Blood pressure range: 30 - 300 mmHg
10. Accuracy of measurement: ± 3 mmHg
11. Range of heartbeat/ pulse rate: 40 - 199 / min
12. Operating conditions - temperature: 5°C ~ 40°C
13. Operating conditions – air humidity: < 80%
14. Storage temperature: -20°C~55°C
15. Conditions for air humidity when storing: < 95%
16. External pressure: atmospheric, 86 kPa ~ 106 kPa
17. Durability of the batteries: about 2 months when using 3 minutes per day
18. Material cuff: Terylene (surface), Nylon (part of sticking), PVC (internal part + tube)

Contents

- 1 **smartLAB[®]pressure** blood pressure monitor
- 1 comfort cuff (size: 22 - 42 cm)
- 1 User manual
- 1 **smartLAB[®]** blood pressure passport
- 4 „AAA“ batteries


General Information:

1. Please relax 5 minutes before measuring.
2. The arm cuff should be placed on top of the elbow (please refer to the indication on the cuff)
3. Do not move or speak when measuring.
4. Please always use the same arm.
5. Please relax 3 to 5 minutes between the measurements so that your blood circulation can start again in your arm.
6. When you do not use the device for one month or more please remove the batteries in order to prevent them from discharging automatically.
7. The device is dedicated for adults with an arm circumference from 22 - 42 cm and not for children or infants.
8. The device works according to its specifications even in case of general cardiac arrhythmia as well as arterial or ventricular cardiac fibrillation.
9. The blood pressure values measured by this device are comparable to the values identified by skilled persons using arm cuff and stethoscope and they range within the required limitations of the American National Standard for electronic or automatic sphygmomanometer.
10. Please avoid strong magnetic interferences such as mobile phone, microwaves etc.

III. Setup & operating functions

Making the device ready for operation

- Please open the cover of the battery case on the bottom side of the device.
- Please insert the „AAA“ batteries according to the required polarity.
- Then close the cover of the battery case.

When the display shows the symbol  for low batteries, please change them all.

Please do not use rechargeable batteries for this device.

When you do not use the device for one month or more please remove the batteries in order to prevent them from discharging automatically.

Measurement device, batteries and arm cuff have to be disposed of according to the local regulation.

Setting time, date and average

1. When the monitor is off, please press " " to enter the year settings

Note:


- During the process of setting, you can press " " button to stop setting at any time.
- If there is no operation during the process of setting, it will turn off within 1 minute.



2. Press “◀ 🔍” button to change the [YEAR]. Each press will increase the numeral by one in a cycling manner. If you hold pressing “◀ 🔍” button, the numeral will increase quickly.


3. Press “⚙️ ▶” to confirm [YEAR], it diverts to [MONTH] setting.

4. Repeat steps 2 and 3 to set [MONTH] and [DAY].

5. Repeat steps 2 and 3 to set [HOUR] and [MINUTE].

6. Repeat steps 2 and 3 to set [UNIT].


7. After the unit is set, the LCD will display "done" and then turn off.

Select the User

1. When the monitor is off, press and hold "Q" button to enter user setting mode. The user ID will blink.

2. Then press "Q" button again, select the user ID between user 1 and user 2.


3. After selecting the suitable user ID, press “▶” button to confirm. The LCD will display “user ID”, and then turn off.


Connecting the cuff to the monitor

Insert the cuff tubing connector into the socket of your **smartLAB[®]pressure** on the back side of the monitor. Make sure that the connector is completely inserted to avoid air leakage during blood pressure measurements.


Applying the cuff

- Place the cuff around a bare arm 1-2cm above the elbow joint.

- While seated, place palm upside in front of you on a flat surface such as a desk or table. Position the air tube in the middle of your arm in line with your middle finger.

- Pull the cuff end through the metal loop (the cuff is packaged like this already), turn it outward (away from your body) and tighten it and close the Velcro fastener.
- The cuff should fit comfortably, yet snugly around your arm. You should be able to insert one finger between your arm and the cuff.

Note:

- Use the same arm for each time you measure.
- Do not move your arm, body, or the monitor and do not move the rubber tube during measurement
- Stay quiet and calm for 5 minutes before blood pressure measurement
- If the cuff becomes dirty, remove it from the monitor and clear it by hand in a mild detergent, then rinse it thoroughly in cold water. Never dry the cuff in clothes dryer or iron it.

Body posture during measurement


Sitting Measurement – Recommended Position

- Sit upright
- Place palm upside in front of you on a flat surface such as a desk or table, with your elbow resting on a chair or table.
- The cuff should be at the same level as your heart.
- Place your arm on the monitor's plastic case as shown in the picture, to aid maintaining a correct and still arm position.


Lying Down Measurement

- Lie on your back.
- Place your arm straight along your side with your palm up-side.
- The cuff should be placed at the same level as your heart
- In case you move your arm and the monitor is not able to measure correctly, the symbol 📶 appears on the display.

Measurement

- After having applied the cuff, seat yourself comfortably and press the  button. All display symbols will be shown.


- After that, the monitor will adjust itself to zero.

- Then, the monitor will start with the measurement by inflating the cuff.

- After the measurement, the air will slowly be released from the cuff. Your blood pressure value and pulse frequency will be shown on the display.


- **smartLAB[®]pressure** will automatically turn off after 1 minute. You can manually shut down the monitor by pressing “ start” after the measurement.
- There are two users in total. Each user has a memory of 60 records.
- You can interrupt the measurement by pressing “ start”. The cuff will then release air and the monitor will turn off.
- If the measurement result is out of the measurement range (SYS: 60mmHg to 230mmHg; or DIA: 40mmHg to 130mmHg; or Pulse: 40-199 pulse/minute), the LCD will display “out”.


Note: The  -symbol indicates that a pulse irregularity has been detected in accordance with the heart beat during the measurement. This is no cause for concern. If the symbol occurs frequently, however, it is advisable to seek medical advice. Please note that the device does not replace a heart examination but is used to detect pulse irregularities at an early stage.

Displaying stored values


- When the monitor is off, please press “ Q” button to show the average value of the latest three records. (Note: If the records are less than 3 groups, the LCD will display the recent record instead. Take user 1 for example.)


- Press “ Q” button or “ ” button to get the records you want.

The order, date and time of the record will be shown alternately.

- If you want to check another user's records, please "start" button to turn off the monitor when the blood pressure monitor is in the memory inquiry mode. Press and hold "Q" button to enter into the selecting user ID mode, press "Q" again to select the user ID between user 1 and user 2, press "G" button to confirm the user ID, then press "Q" button to check the selected user's measurement records.

- You can exit memory mode at any time, by pressing "start". The monitor will then power off.


Note:

The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back on digit (e.g. 2 becomes 3, and so on), and the last record (60) is dropped from the list.

Deleting measured values

If you did not get the correct measurement, you can delete all results by following steps below. (Take User 1 for example.)

- Hold pressing “ Q” button for 3 seconds when the monitor is in the memory recall mode, “User ID+dEL ALL” will show and blink.


- Press “ P” button to confirm deleting and the monitor will display “dEL d0nE” and then turn off.

Note: To exit out of delete mode without deleting any records, press “ start” button before pressing “ P” to confirm any delete commands.


- If there is no record, the following display will show.

Lock the buttons

The unnecessary keys' touch will make the blood pressure monitor turn on and waste electricity. To avoid this, you can press the Lock button to lock the keys if necessary.

- Hold pressing the Lock button until the LCD displays "OFF", this indicates the touch keys have been locked.

- Hold pressing the Lock button until the LCD displays "ON" to unlock the keys. You can press any buttons to use the blood pressure monitor.

IV. Miscellaneous

Maintenance

1. Do not drop this monitor or subject it to strong impact.
2. Avoid high temperature and solarization. Do not immerse the monitor in water as this will cause damage to the monitor. Clean it with a dry cloth in case.
3. If this monitor is stored near freezing, allow it to acclimate to room temperature before use.
4. Do not attempt to disassemble this monitor.
5. If you do not use the monitor for a long time, please remove the batteries.
6. It is recommended to check the monitor performance every 2 years or after repair.
7. If the cuff becomes dirty, remove it from the monitor and clear it by hand in a mild detergent. Do not attempt to clean the reusable cuff with water and never immerse the cuff in water.

Explanation of symbols

Caution! Please read user manual thoroughly!

Type-B Applied Parts


Disposal information

Troubleshooting

Problem	Cause	Remedy
Display will not light up	Batteries are exhausted	Replace with new batteries
	Batteries are inserted incorrectly	Insert the batteries correctly
	AC Adapter is inserted incorrectly	Replace the AC adaptor tightly
LCD is dark / shows „LO“	Battery too weak	Insert new batteries
Display shows „E 01“	Cuff is not secure	Refasten the cuff and measure again
Display shows „E 02“	Cuff too tight	Reapply cuff and start the measurement again
Display shows „E 03“	The pressure of the cuff is excess	Relax for a moment and then measure again
Display shows „E 10“ or „E 11“	The monitor detected motion,talking or the pluse is too poor while measuring.	Movement can affect the measurement.Relax for a moment and then measure again.
Display shows „E 20“	The measurement process does not detect the pulse signal.	Loosen the clothing on the arm and then measure again.
Display shows „E 21“	The treatment of the measurement failed.	Relax for a moment and then measure again.
Display shows „EExx“	A calibration error occurred. (X can be some digital symbol, such as 1, 2,etc., if this similar situation appear, all belong to calibration error.)	Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance.Refer to the warranty for contact information and return instructions.
Display shows „out“	Out of measurement range	Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your physician.


Assessing high blood pressure for adults

The following guidelines for assessing high blood pressure (without regard to age or gender) have been established by the World Health Organization (WHO). Please note that other factors (e.g. diabetes, obesity, smoking, etc.) need to be taken into consideration. Consult your physician for accurate assessment.

Knowledge on Blood Pressure:

Your heart is just like a pump which causes the blood to circulate through all vessels. As a result, the blood pressure is the pressure on the wall of the blood vessel pressed by blood. The output of blood from the heart when it is a systole is called systolic pressure. The blood returning to the heart when it is a diastole is called diastolic pressure.

Blood pressure fluctuates continually — day and night. The highest value usually appears in the daytime and the lowest one usually at midnight. Typically, the value begins to increase at around 3:00AM, and reaches to highest level in the daytime while most people are awake and active.

For this reason it is recommended to measure your blood pressure each day at the same time.

OPERATION & CIRCUMSTANCES	Reason	sytolic pressure		diastolic pressure	
		high	low	high	low
	<i>Drinking, smoking</i>	•		•	
	<i>Excitement or sporting</i>	•		•	
	<i>seeing a doctor, desire to urinate</i>	•		•	
	<i>change of temperature</i>	<i>high</i>	•	•	
		<i>low</i>	•	•	
	<i>After bathing, deep breathing or yawning</i>	•		•	
	<i>Cuff placing</i>	<i>tight</i>		•	•
		<i>loose</i>		•	•

The blood pressure is influenced by a variety of different factors, such as bodily circumstances or a wrong operation of the monitor. Some people even get nervous when seeing a doctor and that will make their blood pressure rise. Single measurements can not indicate the whole blood level, so you need not to be anxious for isolated higher or lower blood pressure results. The most important thing is to master one's blood pressure trend through long-term measurement. Please talk to your doctor for useful judgement of your recorded values.

Average values of normal arterial blood pressure (mmHg):

The following chart shows usual average values. You could obtain similar results by measuring your blood pressure for several days at the same time. (so called „Basic Blood Pressure“)

Age	Male		Female	
	systolic	diastolic	systolic	diastolic
11-15	114	72	109	70
16-20	115	73	110	70
21-25	115	73	110	71
26-30	115	75	112	73
31-35	117	76	114	74
36-40	120	80	116	77
41-45	124	81	122	78
46-50	128	82	128	79
51-55	134	84	134	80
56-60	137	84	139	82
61-65	148	86	145	83

Hypertension:

Various cardiovascular or kidney diseases can cause hypertension. Long-term hypertension can also cause pathological changes in heart, brain and kidney. These changes can lead to hypertension complications which are often diagnosed or treated too late.

Low Blood Pressure:

There are two kinds of minimum pressure:

- If you always feel dizzy, unsuitable in the chest and have difficulties in breathing or occurrences of temporary myopia when standing up, you might have a standing minimum pressure disease. You should go to see a doctor.
- The original minimum pressure disease has no self-symptoms (besides low blood pressure). To get rid of it, you should pay more attention to nutrition and do more physical exercise.


<i>BD Classification</i>	<i>SBP mmHg</i>	<i>DBP mmHg</i>	<i>Color indicator</i>	<i>Color changing area</i>
<i>Optimal</i>	<i>< 120</i>	<i>< 80</i>	<i>GREEN</i>	<i>GREEN</i>
<i>Normal</i>	<i>120-129</i>	<i>80-84</i>	<i>GREEN</i>	<i>GREEN</i>
<i>High-</i>	<i>130-139</i>	<i>85-89</i>	<i>GREEN</i>	<i>GREEN</i>
<i>Stage 1 Hypertension</i>	<i>140-159</i>	<i>90-99</i>	<i>YELLOW</i>	<i>YELLOW</i>
<i>Stage 2 Hypertension</i>	<i>160-179</i>	<i>100-109</i>	<i>ORANGE</i>	<i>RED</i>
<i>Stage 3 Hypertension</i>	<i>≥180</i>	<i>≥110</i>	<i>RED</i>	<i>RED</i>

Note: It is not intended to provide a basis of any type of rush toward emergency conditions/diagnosis based on the color scheme. This color scheme is only meant to discriminate between the different levels of blood pressure.


DPA technology for a more accurate measurement

The **smartLAB[®]pressure** uses a unique measuring technology called DPDA (double pulse detection algorithm) with which the units manage the blood pressure signals twice. The chart left shows the blood pressure signals after managed by DPDA. The result is a more linearized curve, proving the improved measuring accuracy of your **smartLAB[®]pressure**.

Measure-during-inflation

This technology accelerates the measuring process while giving a more comfortable feeling for the users while the cuff is being inflated. It also guarantees a much higher accuracy compared to normal inflating processes.

Regularity and Safety Notice

Safety Information

The signs below might be in the user manual, labelling or other component. They are the requirement of standard and using.

	Caution: Consult accompanying documents		Type BF applied part
CE0537	Symbol for „COMPLIES WITH MDD 93/42/EEC REQUIREMENTS“		DISPOSAL: Do not dispose this product unsorted municipal waste. Collect such waste separately for special treatment is necessary
	Symbol for „MANUFACTURER“		Direct current
	Symbol for „MANUFACTURE DATE“	SN	Symbol for „SERIAL NUMBER“
	Symbol for “Including RF transmitter”		Caution: These notes must be observed to prevent any damage to the device

Complied European Standards List

Risk Management	ISO/EN 14971:2012 Medical devices — Application of risk management to medical devices
Labeling	EN 980:2008 Symbols for use in the labelling of medical devices
User Manual	EN 1041:2008 Information supplied by the manufacturer of medical devices
General Requirements for Safety	EN 60601-1: 2006+A1:2013 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance EN 60601-1-11: 2010 Medical electrical equipment -- Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Electromagnetic Compatibility	EN 60601-1-2:2007 Medical electrical equipment Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests

Performance requirements	EN ISO 81060-1:2012(new) Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type (ISO 81060-1:2007) EN 1060-3:1997+A2:2009 Non-invasive blood pressure Part 3: Supplementary requirements for electromechanical blood pressure measuring system
Clinical investigation	EN 1060-4: 2004 Automatic Blood Pressure Monitor overall system Inter-ventional accuracy of the testing process
Usability	EN 60601-1-6: 2010 Medical electrical equipment -- Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability EN 62366: 2007 Medical devices - Application of usability engineering to medical devices
Software life-cycle processes	EN 62304:2006/AC: 2008 Medical device software - Software life cycle processes

Declaration of Conformity

This product is compliant with the EC Directives: -93/42/EC Annex V

This product is compliant with the EC Directives: -1999/5/EC

EMC Guidance

1. This equipment needs to be installed and put into service in accordance with the information provided in the user manual;
2. Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance $d=3,3m$ away from the equipment.

Note: As indicated in Table above of IEC 60601-1-2:2007 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2 W yields $d=3, 3m$ at an IMMUNITY LEVEL of 3V/m

Safety Statement

This product has been verified to conform to the safety certification in accordance with the specification of EN Standards: EN 60950-1:2006 + A11:2009 + A1:2010 + A12:2011

WEEE note

The WEEE (Waste Electrical and Electronic Equipment) Directive, that came into effect as European law on 13th February 2003, led to a major change in the disposal of electrical equipment.

The WEEE logo on the product and the packaging indicates that it is not allowed to dispose of the product in the ordinary household waste. It is your responsibility to deliver all the disused electrical and electronic equipment to the respective collection points. A separate collection and proper recycling of electronic waste helps dealing economically with natural resources. Furthermore, the recycling of electrical waste is a contribution to the conservation of our environment and thus the health of all people. More information concerning the disposal of electrical and electronic devices, reprocessing and the collection points you can get from the local authorities, waste disposal companies, from retailers and manufacturers of the device.

The primarily purpose of this directive is to prevent electronic waste (WEEE), while encouraging for reuse, recycling and other forms of reprocessing to reduce waste.

RoHS compliance

This product complies with Directive 2011/65/EC of the European Parliament and the Council regarding the limited use of dangerous substances in electrical and electronic equipment (RoHS) and its variations.

Warranty

HMM Diagnostics GmbH products need to fulfill high quality requirements. Because of this reason, HMM Diagnostics GmbH gives a 2-year warranty by purchasing this **smartLAB**[®] product. You can even extend the warranty from 3 to 5 years without extra pay when you register your product. Please register online under following website: www.hmm.info/en/registration.

Wear parts, batteries etc. are excluded from warranty

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