smartLAB[®] sprint nG

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smartLAB[®] sprint nG

Self-Monitoring Blood Glucose System

User Manual

Please read this manual thoroughly before first using this device

HMM Diagnostics GmbH
Friedrichstraße 89
D-69221 Dossenheim, Germany
www.hmm.info

Content

I. Introduction	47
Measuring principle	48
Intended use	48
Important information	49
Note on smartLAB[®] nG „NO CODE“-test strips	51
II. Your smartLAB[®] <i>sprint nG</i>	52
Display & functions	52
Specifications	53
The smartLAB[®] nG test strips	54
Explanation of symbols	55
Set Content	56
III. Setup & functions	57
Insert and changing batteries	57
Setting Time and Date	59
Using smartLAB [®] nG Blood Glucose Test Strips	60
Funktion Check	61
Control solution testing	62
Running a blood glucose test	66
Alternate Site Testing (AST) with the transparent cap	69
Understanding your test results	70
Memory, Storing and displaying stored results	73

IV. Miscellaneous	74
Maintenance	74
Frequently Asked Questions	74
Troubleshooting	77
Limitations of the measurement procedure	80
Regularity and Safety Notice	83
Warranty	83

I. Introduction

Thank you for using the **smartLAB**[®] *sprint nG* Self-Monitoring Blood Glucose System (SMBG). This system was designed to be dependable, easy-to-use, compact, lightweight and portable to help you monitor your blood glucose on a regular basis.

ENGLISH

Please read this manual thoroughly before you begin testing. It provides you and your diabetes care team with important information and step-by-step direction to use the system correctly. Although your **smartLAB**[®] *sprint nG* system is easy to use, you should consult your healthcare professional (this may be your doctor, diabetes nurse educator or pharmacist) for instructions on how to use the system. The only way to obtain accurate results from any system is to make sure you correctly use it. If you have any concerns, please call your authorized distributor in your country during business hours.

Measuring principle

The **smartLAB**[®] *sprint nG* Self-Monitoring Blood Glucose System is designed to allow rapid measurement of blood glucose by using an electrochemical biosensor technology. This system employs a disposable dry reagent strip technology, based on the FAD-binding glucose dehydrogenase.

Each test strip features an electrode containing FAD-binding glucose dehydrogenase. A capillary blood sample is applied to the collecting area of the strip, and is automatically drawn into the reaction zone, where the FAD-binding glucose dehydrogenase catalyzes the glucose dehydrogenation. It will produce gluconolactone. During the reaction, a mediator transfers electrons to the electrode surface and generates a current. The amount of the current is proportional to the amount of glucose present in the blood sample. After 5 seconds, the **smartLAB**[®] *sprint nG* will show glucose concentration on the screen.

Intended use

The **smartLAB**[®] *sprint nG* is intended for the quantitative measurement of glucose in venous whole blood or fresh capillary whole blood from fingertips, palm and forearm. Testing is done outside the body (In vitro diagnostic use). It is indicated for self-testing (over the counter [OTC]) by persons with diabetes, or in clinical settings by healthcare professionals, as an aid to monitor the effectiveness of diabetes control. The meter should be used only with **smartLAB**[®] *nG* blood glucose test strips.

Important information

- The **smartLAB**[®] *sprint nG* blood glucose system is designed and approved for testing venous whole blood or fresh capillary whole blood samples from your fingertips, palm and forearm. The meter is for in vitro diagnostic use ONLY (for testing outside the body). It should not be used to diagnose diabetes.
- The **smartLAB**[®] *sprint nG* blood glucose system can only be used with **smartLAB**[®] *nG* Blood Glucose Test Strips. Other test strips will give inaccurate results.
- Severe impact may cause the meter to malfunction. Do not disassemble the meter as it may damage components inside and cause an incorrect reading. The warranty will be void if the meter has been disassembled.
- Testing is not valid for neonatal blood specimens.
- Incorrect results may occur when performing the test. If you believe you are not feeling well, please contact your healthcare professional immediately.
- Always keep the meter clean and store it in a safe place. Protect from direct sunlight to ensure a longer lifespan.
- The strip slot should be kept free from dirt, dust, blood stains, and water stains.
- Do not store the system and test strips in a car, a bathroom or

a refrigerator. And the system, strips and lancing device should be kept away from children or pets.

- Remove batteries if the meter will not be used for one month or more.
- Store the kits in a dry place with temperature range is 2 to 30°C (35.6 to 86 °F). Keep away from direct sunlight and heat.
- Store your test strips in their original vial only. Do not transfer them to a new vial or any other container.
- Indicate the date you open the vial. D Discard all unused strips immediately after 90 days from open date. The strips are for single use only.
- You should not test critically ill patients with home-use blood glucose meters.
- Warning for potential biohazard: Healthcare professionals using this system on multiple patients should be aware that all products or objects that come in contact with human blood, even after cleaning, should be handled as if capable of transmitting a viral disease.
- Please dispose device according to the local rule of the disposition of electronic device / accessory waste.
- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be

impaired.

-  Do not touch the test strip with wet hands
-  Do not use expired test strips
(see expiration date on the strip vial)
-  Do not bend, cut or twist the strip
-  Consult with your healthcare professional before testing
on your palm or forearm.
-  Altitude up to 3,048 meters above sea level has no effect
on readings.

Health-Related Information

- If you are experiencing dehydration, frequent urination, low blood pressure, shock or hyperosmolar hyperglycemic nonketotic coma (HHNKC), you may get a test result that is lower than what your blood glucose really is. If you think you are dehydrated, call your healthcare professional right away.
- If you have followed the steps in the user's manual, but still have symptoms that do not seem to match your test results, or if you have questions, please contact your healthcare professional.
- Please read your test strip instructions carefully for additional health-related information

Note on smartLAB[®] nG „NO CODE“-test strips

The **smartLAB[®] nG** blood glucose test strips compatible with **smartLAB[®] nG** glucose meters need not to be coded, any more.

II. Your smartLAB® *sprint* nG

Display & functions

Folder Symbol		Meaning
88-88	Date	Show the moment of month and day.
88:88	Time	Show the moment of time.
888	Test result	Show the test result.
888	Record/Average display	Tracking all the test results or average of /14/30/90 days.
CS	Control Solution	Control Solution test and the result is not included in the memory.
A	AM	Indicate before noon
P	PM	Indicate after noon
	Batteries status	When the battery symbol appears, prepare new batteries for installation.
	AC = ante cibos	Indicate before meal
	PC = post cibos	Indicate after meal
	Insert strip	Insert the test strip to test blood glucose.
	Thermometer	When the thermometer icon appears, the temperature is too high or too low.
	Error	The meter system is in fault.
	Memory	Show the memorized results
	Apply blood	Apply blood to the test strip
mg/dL	Glucose unit	The testing unit of blood glucose is mg/dL.
mmol/L	Glucose unit	The testing unit of blood glucose is mmol/L.

Specifications

1. *Type:* smartLAB® *sprint* nG
2. *Measuring range:* 20~630 mg/dL (1.1~35.0 mmol/L)
3. *Measuring duration:* 5 seconds
4. *Time mode:* AM/PM - 24 hours
5. *Memory:* 480 values
6. *Operating temperature:* 10°C~40°C (50°F~104°F)
7. *Storage Temperature:* -20°C ~ 50°C (-4°F~122°F)
8. *Relative humidity:* RH ≤90%

9. *Blood sample:* ≥ 0,6 µL capillary fresh blood or venous whole blood from fingertips, palm, or forearm
10. *Calibration:* Plasma equivalent
11. *Hematocrit (Hct):* 30-55%
12. *Battery type:* 2 1,5 V AAA Lithium batteries
13. *Battery life:* over 2000 measurements
14. *Display-size:* 43 x 35 mm
15. *Meter dimensions:* 98 x 57 x 22 mm
13. *Weight:* 53,2 g (without batteries)


Front

Test Strip Slot-
Insert test strip here.

Display-
Shows results, messages, and results stored in memory.


Left Button-
Press to enter memory, adjust setting, and scroll through results.

Right Button-
Press to turn the meter on, confirm setting, and scroll through results.


Side

Strip Ejector-
Push the Strip Ejector to remove the strip.


Back

Battery Door-
Flip open the battery door by pushing the tab in the direction of the arrow and pulling the door up.



The smartLAB[®] nG test strips

Note: You can find the test strip expiration date on the vial label next to the ⌚ symbol.

Explanation of symbols


-  *Caution!*
-  *Lot number*
-  *Expiration date*
-  *For in vitro diagnostic use only*
-  *1.5V(AAA) x 2 batteries only*
-  *Single use only*
-  *Storage temperature*
-  *Manufacturer*
-  *Before use, read manual first*
- REF** *Catalogue number*
- SN** *Serial number*
-  *Control solution*
-  *Opening date of test strip vial*
-  *Sufficient for*
-  *Blood glucose test result in mg/dL*
-  *Blood glucose test result in mmol/L*
-  *Green Dot / Duales System Deutschland GmbH (DSD)*
-  *This product meets the requirements of Directive 98/79/CE for in vitro diagnostic medical devices.*
-  *5-seconds result*
-  *Keep away from sunlight*
-  *0.6 µL blood volume*
-  *Human-factor flow*
-  *Large LCD screen*
-  *Discard 6 months (180 days) after opening*

Set contents

- 1 **smartLAB**[®] *sprint nG* blood glucose meter
- 1 User manual
- 1 carrying bag
- 2 AAA 1,5V Alkali batteries
- 1 **smartLAB**[®] lancing device
- 1 **smartLAB**[®] Check Strip
- 1 **smartLAB**[®] *nG* Quick Start instruction
- 10 **smartLAB**[®] lancets
- 10 **smartLAB**[®] *nG* blood glucose test strips
- 1 **smartLAB**[®] *nG* blood glucose test strip instruction

III. Setup & functions

Insert and changing batteries

1. Open the battery door on the back of the meter by pushing the tab in the direction of the arrow and pulling the door up.
2. Insert two batteries and the meter beeps.
3. Put the battery door back in place and snap it closed. The meter turns on automatically.

Note:

1. Remove the batteries when you will not be using the meter for one month or longer.
2. Values won't be deleted when the batteries are replaced.


Meter, batteries, lancets, test strips etc. must be disposed of according to local regulations at the end of their usage.

Low Battery

The meter uses two alkaline 1.5V (AAA) batteries. Batteries will normally last for more than 2000 tests. Other types of 1.5V (AAA) batteries are also acceptable, but the capacity of test times may differ. Insert the batteries when you first use the meter or replace with new batteries when the “LP” (low power) message and the low battery symbol appear on the display.

The meter will not turn on the first time batteries are inserted. Please press and hold right button or insert the test strip to turn your meter on.

The meter will turn off automatically. Or you can press and hold right button to turn your meter off.

Note:

1. The meter won't delete earlier records after you replace batteries.
2. You should reset the time and date again after you replace the batteries.
3. 1.5V (AAA) x 2 batteries are available at most stores. You may take the old batteries with you for replacement.
4. Remove batteries when you will not be using the meter for one month or more.

Setting Time and Date

Setting the current time and date in your meter is important if you use the meter memory.


1. Press the right button to turn the meter on.

2. The display shows the last 2-digit of the year that flashes at the top of display as well. Press the left button to adjust the year and press the right button to confirm the setting.

3. Repeat step 2 to set date and time. The flashing field is the one you are currently setting. (A=AM, P=PM)

Note: AM means the time from 00:00 o'clock until 12:00 o'clock. PM means the time from 12:00 o'clock until 24:00 o'clock.

Using smartLAB[®] nG Blood Glucose Test Strips

- Use only with **smartLAB[®] nG** Blood Glucose Systems.
- Run a control solution test every time you open a new box of test strips (See Chapter „Control Solution Testing.“)
- Keep the test strips in their original bottle.
- After you take a test strip out of the bottle, tightly close the bottle immediately. This keeps the test strips dry.
- Use the test strip within three minutes after taking it out of the bottle.
- The strip is for single use only. Do not reuse it.
- Record the date you open the test strip bottle. Be sure to check the expiration date on the test strip bottle. The test strip is good for 6 months from the date the bottle is opened or until the expiration date on the bottle, whichever comes first.
- Store the test strip bottle and your system in a cool dry place.
- Store the test strips between 2°C -30°C (36°F - 86°F). Do not freeze.
- Do not apply blood or control solution to the test strip before you insert it into the system.
- Do not touch the test strip with wet hands. Do not bend, cut, or twist the test strips.

Function check

You can carry out a function check with the Check Strip which is included in your system set. You do not have to do the function check before every measurement. It helps you to control if your **smartLAB[®]nG** system works properly, from time to time, though.

1. Insert check strip into the meter:

Make sure your check strip is inserted with writing to top. The meter will turn on automatically and the display will show „CHK“.

2. Check message in LCD display:

The function check will take approx. 3 seconds. The check can produce the following messages on the display:

- „OK“ - your meter works properly
- „FAL“- change batteries and check again. If the message still appears, please refer to your local supplier.

3. Remove check strip:

The meter will turn off automatically after removing the test strip. Please keep the check strip in the soft case of the meter.

Control Solution Testing

Running a control test lets you know that your meter and test strips are working properly to give reliable results. You should run a control test when:

- You use the **smartLAB[®] *sprint nG*** Blood Glucose System for the first time.
- You open a new vial of test strips.
- You think the meter or test strips might be working incorrectly.
- You drop the meter.
- You have repeated a test and the test results are still lower or higher than expected.
- You are practicing the test procedure.

Note: Professional users are instructed to follow federal, state, and local guidelines concerning QC practices.


About the **smartLAB[®] *nG*** control solution

- **smartLAB[®] *nG*** control solution is not intended for human consumption or injection.
- Use only with **smartLAB[®] *nG*** test strips.
- Write the date you first opened the bottle on the bottle label.
- The **smartLAB[®] *nG*** control solution is durable for three months from the date the vial is opened or until the „Expiration Date“ on the bottle, whichever comes first.
- Do not use **smartLAB[®] *nG*** control solution that is past the „Expiration Date“

- The **smartLAB[®]nG** control solution can stain clothing. If you spill it, wash your clothes with soap and water.
- Close the bottle tightly after use.
- Left over control solution should not be added back into the control bottle.
- Store the bottle of **smartLAB[®]nG** control solution at room temperature, below 30°C (86°F). Do not freeze nor refrigerate.
- If you would like to purchase **smartLAB[®]nG** control solutions, please contact your local dealer.

Running a smartLAB[®]nG control solution test

Please make sure you have all necessary items at hand when running a control solution test. You need the **smartLAB[®]sprint nG** meter, a **smartLAB[®]nG** test strip, and **smartLAB[®]nG** control solution.

1. Put a test strip into the meter as the direction of the arrow on the test strip. The meter turns on automatically.

2. Press the left button to set  (AC),  (PC) or CS Mode, then press right button to confirm the setting. Select CD Mode here.

3. The icons  (insert strip) and  (apply blood) show themselves.


4. Place the meter on a flat surface, like a table.

5. Remove the control solution bottle cap and wipe the tip of the bottle with a tissue.

6. Squeeze the bottle until a tiny drop forms at the tip of the bottle.

7. Touch the drop to the Blood collection area at the tip of the test strip. Do not put control solution on top of the test strip.

8. The meter starts to count down from 5 seconds and then shows the test result.

9. Do not remove the test strip yet. Check if the reading falls within the range printed on the test strip bottle.

10. Push the Strip Ejector to remove the test strip.

11. Throw the test strip away after you have compared the reading to the range printed on the test strip bottle.

Note: Be aware that there are Normal and High ranges listed on your test strip bottle. The Normal and High refer to the two different control solutions that are available. Please be sure that you are looking at the correct range.

4. Understanding control test results

The label on your test strip vial shows the acceptable ranges for the **smartLAB[®]nG** control solutions. The result you get should be inside this range. Make sure you compare the result to the correct level of control. When the control result is inside the range on the test strip vial, your test strips and your meter are working properly. If your control result is not inside the acceptable range (printed on your test strip vial), here are some things you can do to solve the problem:

Problem	Solution
<i>Was the test strip exposed to open air for a long period of time?</i>	<i>Repeat the control test with properly stored strips.</i>
<i>Was the test strip vial capped tightly? Or was test strip cap left open?</i>	<i>If the cap was not tight, or the bottle was left uncapped, open a new bottle of test strips. Do not reuse the strips from the affected bottle.</i>
<i>Does the meter work properly?</i>	<i>You can use the check strip to verify the system's functions.</i>
<i>Is the control solution expired or contaminated?</i>	<i>Replace with new control solution to check the performance of your glucose system.</i>
<i>Were test strips and control solution stored in a cool and dry place?</i>	<i>Repeat the control test with properly stored strips or control solutions.</i>
<i>Did you follow the testing steps properly?</i>	<i>Read Chapter "Control Solution Testing" again and retest.</i>

Note:

1. If "[5]" is not set, then the control test values will be counted into 7/14/30/90 days average mode.
2. If "[5]" is set, the control solution values will not be counted into averages, but will be stored in the memory for review.


Running a blood glucose test

1. Preparations

Before running a blood glucose test, make yourself acquainted with the test strips and the lancing device. Keep all required materials ready for use: your **smartLAB[®]sprint nG** meter, the **smartLAB[®]nG** test strips and the **smartLAB[®]** lancing device including lancets. Wash your hands thoroughly with warm water before collecting the blood sample. Rinse your hands thoroughly.

2. Collecting a blood sample

When the blood drop flashes on the display, gently massage the part of your finger tip which you obtain a drop of blood from using the lancing device. Place the lancing device against the pad of your finger. Press the trigger button to activate the lancing device.


3. Blood glucose measurement


1. Put a test strip into the meter in the direction of the arrow. The meter turns on automatically.

2. Press the left button to set 🍏 (AC) or 🍏 (PC) Mode, then press the right button to confirm the setting.

3. The icons 📄 (insert strip) and 💧 (apply blood) show themselves. When the blood drop flashes on the display, obtain a drop of blood from your finger.

4. Hold the lancing device firmly against the side of your finger. Press the trigger button. Then remove the lancing device from your finger.
5. Gently squeeze and/or massage your fingertip until a round drop of blood on your fingertip appears. Wipe away the first drop with a tissue and use the second drop.
6. Touch the drop to the tip of the transparent window of the test strip. Do not put blood on top of the strip. Be sure to get enough blood in the strip's confirmation window. Otherwise, an inaccurate reading may yield.

1. The system starts to count down from 5 seconds and then shows the test result.
2. Push the Strip Ejector to remove the test strip.
3. Throw it away after you have compared the reading to the range printed on the test strip bottle.

Note: To avoid infections, never let anyone else use your lancing device and lancets. Sterile lancets are for single use only.

 Please make sure to apply blood when the blood drop appears on the display.



 Do not reuse lancets.

 Dispose lancets and test strips according to your local regulations.

4. Discarding used test strips and lancets

Write down the measured result in your personal diabetic log book and take the test strip off the strip slot. The meter turns off. If you don't take off the test strip, the meter will turn off automatically after 5 minutes.

In order to discard lancets from the **smartLAB**[®] lancing device, remove the end cap and carefully pull the used lancet out of the holder. Dispose the lancet according to local regulations, in order to avoid hurting other people.

Alternate Site Testing (AST) with the transparent cap

You can also gain the blood sample from other parts of the body than the finger tip. If you want to use Alternate Site Testing, please use the transparent cap for your lancing device. Carry out the following steps:

- Gently massage the desired punctuation site on your arm or hand for a few seconds. This will improve the blood circulation.
- Push the lancing device with the transparent cap against the desired blood collection site and push the trigger button to gain a sample.
- Exercise a constant pressure until you see through the transparent cap, that enough blood is gained. After that you can carry out the blood glucose test, as usual.

Alternate Blood Collection Sites

Understanding your test results

The **smartLAB[®]nG** blood glucose test strips are plasma referenced and calibrated for easier comparison to lab results. The unit of blood glucose test results displayed on the screen is either mg/dL or mmol/L, depending on which unit of measurement you have selected. The mmol/L results will always include a decimal point; mg/dL results do not include a decimal point.

Expected results for non-diabetic adults:

The normal fasting glucose range is 70 to 100 mg/dL (3.9 to 6.1 mmol/L). Two hours after meals, normal glucose values should be less than 120 mg/dL (6.7 mmol/L).

Expected results for diabetic adults:

The normal fasting blood glucose range for an adult with diabetes is 70-130 mg/dL (3.9-7.2 mmol/L)*. Two hours after meals, the blood glucose range for an adult with diabetes is less than 180 mg/dL (10.0 mmol/L). For further queries about diabetes: please consult your healthcare professional for the blood glucose range appropriate for you.

Unusual test results:

If your blood glucose result doesn't match the way you feel, follow these steps, and then repeat the test:

1. Check if the strips are within the expiration date.
2. Be sure that the drop of blood completely filled the reaction zone of the test strip

*Reference: American Diabetes Association. Standards of medical care in diabetes-Table 10. Diabetes care. 2011; Vol. 34, Suppl 1, S21.

3. Do a function control test with the **smartLAB**[®] check strip
4. *Optional:* Check meter and test strip performance with the **smartLAB**[®] *nG* control solution. When test results are still questionable or inconsistent, consult your healthcare professional before making any changes to your diabetes medication program.

5. Repeat the Measurement:

If the values of the repeated measurements are however questionable and do not match the way you feel, please contact your health care professional before you change your medication.

Note:

1. Extremely high humidity may affect the test results. A relative humidity greater than 90% may cause inaccurate results.
2. A red blood cell count (Hematocrit) below 30% may cause higher results. Hematocrit above 55% may cause lower results.
3. Some studies have shown that electromagnetic fields may affect results. Do not test near an operating microwave oven.

Symptoms of high or low blood glucose:

Being aware of the symptoms of high or low blood glucose can help you understand your test results and decide what to do if they seem unusual. A high blood glucose would be greater than 240 mg/dL (13.33 mmol/L). A low blood glucose would be below 60 mg/dL (3.33 mmol/L). Here are the most common symptoms:

High blood glucose (hyperglycemia): fatigue, increased appetite or thirst, frequent urination, blurred vision, headache, general aching, or vomiting.

Low blood glucose (hypoglycemia): sweating, trembling, blurred vision, rapid heartbeat, tingling, or numbness around mouth or fingertips.

If you are experiencing any of these symptoms, test your blood glucose. If your blood glucose result is displayed greater than 240 mg/dL (13.33 mmol/L) or below 60 mg/dL (3.33 mmol/L) and you have symptoms of low or high blood glucose, contact your doctor immediately. If your blood glucose result does not match how you feel, follow the steps under „Unusual Test Results.“

Comparing your meter's result to a lab result:

A common question is how the blood glucose results on your meter compare to the lab results. Your blood glucose can change quickly, especially after eating, taking medication, or exercising. If you test yourself in the morning, then go to the doctor's office for a blood glucose test. The results will probably not match, even if you are fasting. This is typically not a problem with your meter, it just means that time has elapsed and your blood glucose has changed.

If you want to compare your meter result to the lab result, you must be fasting. Bring your meter to the doctor's office, and test yourself by fingertips within five minutes of having blood drawn from your arm by a healthcare professional. Keep in mind that the lab could use different technology than **smartLAB®** *sprint nG* blood glucose meter, and that blood glucose meters for self testing generally read somewhat lower or higher than the lab result.

For accuracy and precision data and for important information on limitations, see the instructions that come with your test strips.

Memory, Storing and displaying stored results


Your meter stores a maximum of 480 test results with the time and date of the test. You can review them at any time. When the memory is full, the oldest result is dropped as the latest is added, so it is very important to have the correct time and date set in the meter.

Note:

1. Do not change your therapy based on one individual result in memory.
2. The memory is not lost when you replace the battery. You do need to check that the time and date are still correct for future readings. See Chapter „Setting the time and date“.
3. Once 480 results are in memory, adding a new result causes the oldest one to be deleted.

Viewing Test Results

The meter provides 7, 14, 30 and 90 days averaging to help track your blood glucose result soundlessly.

1. Press the right button to turn on the system.
2. Press left button to view the latest test result.
3. All results are shown in turn from records of M480 to M001.
4. Press the right button to view averages of 30day-AC/PC and 7/14/30/90 days.

IV. Miscellaneous

Maintenance

Your **smartLAB[®] *sprint nG*** blood glucose meter does not require any special cleaning. Just keep the meter free of dirt, dust, blood- and water stains. Following these guidelines carefully will help you getting the best performance possible: Gently wipe the meter's surface with a soft cloth Do not get any moisture in the test strip slot.

Please do not spray any cleaning solution directly onto the meter.

Do not pour liquid into the meter or put the meter under water or liquid.

To clean the lancing device, use a mild dishwashing liquid and a soft cloth. **DO NOT** place the entire device under water.

Please make sure that no liquids enter the meter's apertures.

Frequently Asked Questions


1. The device won't turn on. Please check:
 - Are the batteries inserted correctly?
 - Are the batteries empty and need to be replaced?

2. No symbol after inserting the test strip. Please check:
 - Did you insert the test strip correctly? The strip needs to be inserted arrow up in the direction of the arrow.
3. No symbol after inserting the check strip. Please check:
 - Did you insert the dark-blue check strip correctly? The writing must be to the top.
4. No measurement after blood sample has been sucked in. Please check:
 - Is the reaction zone filled completely (red)? Repeat the glucose test with a new strip.
 - Did you insert the test strip correctly (see 2.)?
 - Does the display read „ERR“? The inserted test strip was already used. Please use a new test strip.
5. Blood glucose result too high / low. Please check:
 - Was the test strip exposed to air for a long time? Take out the strips from the vial directly before the measurement, since humidity and other environmental influences may affect the test result.
 - Was the test strip vial tightly closed? Humidity and other environmental influences may affect the test results!
 - Alternatively, check the function of your device with the check strip. If after 3 seconds the message “OK” appears on the display, the device is working properly.
If this is not the case, please contact your local supplier.

6. Hematocrit: Test strip results are not significantly affected by hematocrits in range of 30% to 55%. Hematocrit levels less than 30% may cause incorrect high readings and hematocrit levels greater than 55% may cause incorrect low readings. If you do not know your hematocrit level, consult your healthcare professional.
7. Patients undergoing oxygen therapy may have inaccurate results.
8. Altitude up to 3050 meters above sea level has no effect on readings.
9. Test results may be false if the patient is severely dehydrated or severely hypertensive, in shock, or in hypoglycemic-hyperosmolar state (with or without ketosis). Critically ill patients should not be tested with home-use blood glucose meter.
10. Increased cholesterol and triglyceride levels may interfere with the way light is reflected producing erroneous meter results.

Troubleshooting

Note: If you are not sure how to react on error codes, please contact your local supplier.

Battery empty

Display: „LP“ & „Battery symbol “

Solution: Replace with new batteries.

Test strip already used or wet

Display: „Err“ & „Insert“

Solution: Use a new test strip.

System error

Display: „001 Error“

Solution: Replace batteries first. If the error still occurs, please contact your local supplier.

System Error

Display: „005“

Solution: Replace the batteries first. If Err appears again, please contact your local dealer.

Test result is higher than 630 mg/dL (35.0 mmol/L)

Display: „HI“

Solution: Test again. If the result is still the same, try a control solution test and if the control solution test falls within the correct range, please call your healthcare professional immediately.

Test result is under 20 mg/dL (1.1 mmol/L)

Display: „LO“


Solution: Test again. If the result is still the same, try a control solution test and if the control solution test falls within the correct range, please call your healthcare professional immediately.

Err message. Not enough sample was drawn into the test strip for measurement.

Display: „Err & Blood Drop“

Solution: Discard the test strip and repeat the test.


Temperature too high

Display: „Ht“ & „Thermometer symbol“

The operating temperature is too high (above the required temperature range from 10°C - 40°C (50°F - 104°F)). This alerts users that an incorrect result may occur if the test continues.

Solution: Relocate the meter to a location with temperature between 10°C - 40°C (50°F - 104°F).

Temperature too low

Display: „Lt“ & „Thermometer symbol“

The operating temperature is too low (below the required temperature range from 10°C - 40°C (50°F - 104°F)). This alerts users that an incorrect result may occur if the test continues.

Solution: Relocate the meter to a location with temperature between 10°C - 40°C (50°F - 104°F).

Limitations of the measurement procedure

1. DO NOT use serum or plasma sample.
2. DO NOT use neonate blood sample.
3. DO NOT use anticoagulant NaF or potassium oxalate for venous sample preparation.
4. Extreme humidity may affect the results. A relative humidity greater than 90% may cause incorrect results.
5. The system is designed to be used at temperatures between 10°C and 40°C (50°F and 104°F). Outside this range, the system may yield erroneous results.
6. DO NOT reuse the test strips. The test strips are intended for single use only.
7. DO NOT use iodoacetic acid, fluoride or sodium fluoride/oxalate as a preservative for blood specimens.
8. Hematocrit:
Test strip results are not significantly affected by hematocrits in range of 30% to 55%. Hematocrit levels less than 30% may cause incorrect high readings and hematocrit levels greater than 55% may cause incorrect low readings. If you do not know your hematocrit level, consult your healthcare professional.
9. Critically ill patients should not be tested with **smartLAB[®]nG** Blood Glucose Systems .
10. DO NOT use during xylose absorption testing. Xylose in the blood will interfere Self-Monitoring Blood Glucose System.

11. Patients undergoing oxygen therapy may have inaccurate results.
12. Altitude up to 3050 meters above sea level has no effect on readings.
13. Test results may be false if the patient is severely dehydrated or severely hypertensive, in shock, or in hypoglycemic-hyperosmolar state (with or without ketosis). Critically ill patients should not be tested with home-use blood glucose meter.
14. Elevated cholesterol and triglyceride levels may interfere with the way light is reflected producing erroneous meter results.
15. Recent studies have shown that EMI can cause electronic medical device performance degradation and could lead to inappropriate therapy.
16. Lipemic samples: Cholesterol level up to 500 mg/dL (12.92 mmol/L) and triglycerides up to 3,000 mg/dL (33.6mmol/L) do not affect the results. Grossly lipemic patient samples have not been tested and are not recommended for testing with **smartLAB®** *sprint nG* Blood Glucose Meter.

17. Interfering Substances depend on the concentration. The below substances up to the test concentration will not affect the test results:

Concentrations of the interference tested	Bias		Glucose Level	80 mg/dL (4.4 mmol/L)	250 mg/dL (13.9 mmol/L)	500 mg/dL (27.8 mmol/L)
Ascorbic Acid	4 mg/dL	(0.26 mmol/L)		10.89%	-1.76%	4.55%
Ibuprofen	50 mg/dL	(2.43 mmol/L)		3.10%	2.88%	4.62%
L-Dopa	1.8 mg/dL	(0.09 mmol/L)		10.59%	7.91%	4.90%
Sodium Salicylate	50 mg/dL	(3.12 mmol/L)		-2.59%	9.42%	-0.84%
Tetracycline	1.5 mg/dL	(0.03 mmol/L)		-5.32%	3.81%	3.20%
Tolbutamide	100 mg/dL	(3.70 mmol/L)		-2.60%	12.30%	0.89%
Bilirubin-unconjugated	2.4 mg/dL	(0.04 mmol/L)		-2.52%	4.05%	-0.23%
Uric acid	8 mg/dL	(0.48 mmol/L)		2.71%	9.55%	-1.75%
Xylose	4 mg/dL	(0.27 mmol/L)		-5.12%	-1.64%	-4.44%

Regularity and Safety Notice

WEEE note

The WEEE (Waste Electrical and Electronic Equipment) Directive, that came into effect as European law on 13th February 2003, led to a major change in the disposal of electrical equipment.

The primary purpose of this directive is to prevent electronic waste (WEEE), while encouraging for reuse, recycling and other forms of reprocessing to reduce waste.

The WEEE logo on the product and the packaging indicates that it is not allowed to dispose of the product in the ordinary household waste. It is your responsibility to deliver all the disused electrical and electronic equipment to the respective collection points. A separate collection and proper recycling of electronic waste helps dealing economically with natural resources. Furthermore, the recycling of electrical waste is a contribution to the conservation of our environment and thus the health of all people. More information concerning the disposal of electrical and electronic devices, reprocessing and the collection points you can get from the local authorities, waste disposal companies, from retailers and manufacturers of the device.

Warranty

HMM Diagnostics GmbH products need to fulfill high quality requirements.

Because of this reason, HMM Diagnostics GmbH gives a 2-year warranty by purchasing this **smartLAB**[®] product. You can even extend the warranty from 2 to 5 years without extra pay when you register your product. Please register online under following website: www.hmm.info/en/registration.

Wear parts, batteries etc. are excluded from warranty.

Manufacturer:

HMM Diagnostics GmbH
Friedrichstr. 89
D-69221 Dossenheim, Germany

mail: info@hmm.info
www.hmm.info

Weitere Informationen zu den smartLAB® Produkten /
More information on our smartLAB® products:

www.smartlab.org